

# Silicone Manual Resuscitators

**CE**  
2195



**Infant (250 ml)**  
with size 00 mask



**Child (500 ml)**  
with size 01 mask



**Adult (1600 ml)**  
with size 5 mask

## Features

- Rugged; 100% autoclavable; reusable
- Non-rebreathing valve with 40 cms of H<sub>2</sub>O pressure release (except adult)
- 360° swiveling patient connector
- PVC oxygen reservoir bag
- 1.5 mtr PVC oxygen tubing
- With carry bag

## Face Masks



## Oxygen Hood

- Single piece, round shape, unbreakable, autoclavable polycarbonate material
- Trauma-free silicone neck adjustable flap
- Bilateral oxygen nozzle (prevents direct flow of cold oxygen to patient's head)
- Width : 8.5"; Height : 7.5" (Medium), 6.5" (Small)



**Bird Meditech**

ISO 13485:2003 CERTIFIED COMPANY

109, 110, Nidhi Industrial Estate, Shiv Shankar Industrial Complex No. 2,  
NH No. 8, Village Waliv, Vasai (East) 401208, Dist. Thane, Maharashtra, India.  
Tel. : +91 250 321 2729, E mail : birdmeditech@sify.com  
[www.birdmeditech.net](http://www.birdmeditech.net)