




ULTRASONIC SCALPEL

Adoptive Tissue
Technology for
Better Surgical
Precision &
Performance



Ultrasonic Generator using Frequency 50,500 Hz

Reliably seal vessels up to 5mm




Burst Pressure: EG US

Thermal Tissue Damage: <0.3mm



EG US 36




**Superior ergonomics,
Intuitive design Integrated
clicker while closed tissue**

**Precision Coagulation
Minimal Lateral
Thermal Tissue**

**Multifunctional-Grasp,
Dissect, Cut, Coagulate
with a single instrument**

- Tissue does not slip with the tighter gripping
- Improve performance in thicknesses
- Reduce the use of SULK, save the cost
- Generates less stress on tissue during compression and gripping
- Wider proximal-to-distal jaw opening for easy tissue positioning and manipulation.
- 90 degrees turn around the joint, easy to reach the complex surgical site, help to complete the difficult operation



EG US 36R5mm 36cm Curve

EG US 36R5mm 23cm Curve

EG US 14R5mm 14cm Curve



**Marketed by :
Medical Systems & Services**

101, Express Building, H-Block Market,

Ashok Vihar, Delhi-110052

Landline : 011-41080676 / 47091476

Email : mediservesales@gmail.com

Web : www.mediservegroup.com

CE