

Accuracy at it's Best

Digital Weighing Scales



Digital Baby Weighing Scale

20 Kg with ± 5 gm Accuracy

05 Kg with ± 1 gm Accuracy

- Durable, Moulded Plastic Baby Tray
- Feather Touch Keys
- Zero Setting Facility
- Auto Hold Facility for Accurate Weighing of Active Babies
- Kg to Pound Conversion Facility
- User Friendly Calibration Facility
- Operating Voltage : 230V - 50Hz AC



Optional
Battery Backup