

# NASAL CANNULA

## Order Information

With HUMIDOFLO, patients can eat, drink and communicate during therapy. This allows comfortable, effective delivery of humidified, continuous high flow gas, blended with Oxygen up to 100%, creating an ideal solution for a wide range of respiratory failure patients.


Ultra-soft nasal prongs.

Button clip for extra security.

Easy glide neck strap.

Acute hypoxemic respiratory failure

- Post-extubation period
- Pre-intubation
- Emergency department
- Bronchoscopy and other invasive procedures
- Palliative care
- Acute heart failure
- Chronic airway disease



### INFANT NASAL CANNULA

REF	HFN-1001(Pink)	HFN-1002(Green)
Cannula Size	2.6 mm OD	3.5mm OD
Flow Range	10Lpm Max	15Lpm Max

### NASAL CANNULA

REF	HFN-0001	HFN-0002	HFN-0003
Model	Pediatric	Adult/Child	Adult
Color	White	Light Purple	Purple
Size	#S	#M	#L
Weight	40g	40g	40g
Cannula Size	4mm OD	5mm OD	6mm OD
Flow Range	35Lpm	45Lpm	60Lpm
Package	400 sets(8boxes x 50sets) L55 x W55 x H55 cm		

### TRACHEOSTOMY INTERFACE

REF	HFN-0006 (Purple)	HFN-0007(Violet)
Size	For Infant (flow rate 5~15 L/min Prong size 2)	For Adult (flow rate 15~50 L/min) Prong size 5



GREAT GROUP MEDICAL CO., LTD.

No.168, Xingong 2nd Rd., Tianzhong Township,  
Changhua County, Taiwan

Tel:886-4-875-8161 Fax:886-4-875-8161



CE 0434



info@greatgroup.com.tw

www.greatgroup.com.tw